

Abstract

RetroNectin[®] (a recombinant human fibronectin fragment CH-296) significantly enhances the efficiency of viral-mediated gene transfer, and is widely used in various clinical applications involving gene-engineered cells, including T cells and hematopoietic cells.

However, there are several challenges in manufacturing clinical-grade gene-engineered T-cell products that must be addressed in order to standardize manufacturing methods and achieve consistent results. To simplify and streamline the gene transduction process, it is feasible to incorporate RetroNectin as a transduction enhancer directly into the liquid. This modification would reduce the number of steps involved and can be easily implemented in automated cell processing systems.

In this study, we evaluated the efficiency of enhancing gene transduction by directly adding RetroNectin into the liquid phase during the transduction process of T cells using lentiviral vectors. The results showed that the direct addition of RetroNectin and lentiviruses significantly improved transduction efficiency in T cells, with a 1.5 to more than 3-fold increase compared to viral infection without RetroNectin. Importantly, the addition of RetroNectin did not have any adverse effects on cell viability, cell growth, or cell phenotype.

The Structure of RetroNectin[®] and the mechanism of RetroNectin-mediated transduction

RetroNectin[®]

Recombinant human fibronectin fragment CH-296

functional domains:

- Heparin-binding domain(H-domain) \Rightarrow bind to viral particles
- Cell-binding domain (C-domain) ⇒bind to VLA-5*
- · CS-1 sequence \Rightarrow bind to VLA-4*



Enhance lentiviral- and retroviral-mediated gene transduction by facilitating the colocalization of target cells and viral particles.









The standard protocols for RetroNectin mediated transduction

The Direct Addition of RetroNectin[®] Solutions in the Liquid Phase Improves the Efficiency of Lentiviral Gene Transduction.

Yu Okubo, Yoko Kudo, and Sachiko Okamoto,

T cell transduction using various protocols

Transduction : 48-well plate (2x10⁵ cells/0.3mL/well)



Technology Development Center, TAKARA BIO INC., Shiga, Japan

Exp.2: Concentration of RetroNectin & Pre-incubation

Day 0 CD4/8 selection & stimulation (aCD3/CD28 beads)

Day 1 LV Td (VSVG-LV-ZsGreen1) 48 well plate/4x10⁵ cells/mL/0.5 mL/well w/ or w/o pre-incubation of LV&RN 37°C 1h

Day 2 Expansion culture 1:10 Day 5/6 FC Analysis

Pre-incubation slightly increased the Td efficiencies in some conditions, but not significant.

A correlation was observed between the enhancing effect and RetroNectin concentration.

Conclusion

Direct addition of RetroNectin solution can improve the Lentiviral Td efficiencies in T cells

- 1.5 to more than 3-fold increase compared to viral infection without RetroNectin
- Increasing the RN concentration leads to a greater enhancement effect
- No adverse effects on cell viability, cell growth, or cell phenotype

New RetroNectin protocol is simple & can be easily integrated into automated cell processing systems.

> COI: Y.O., Y.K., and S.O are employees of TakaraBio Inc.